

# Over 50 years

The longest established cleanroom company in the UK



## BASSAIRE

IS CLEAN AIR

## Isolators and gloveboxes

### Custom isolators

Bassaire design and produce isolators and gloveboxes to suit a variety of unique situations where potent or toxic hazards need to be addressed.

High integrity isolator designs are suitable for low molecular contamination down to 1-2 ppm oxygen 3-5 ppm moisture.

Materials employed in the construction of these include fully welded 4mm thick stainless steel or lightweight PVC.

### Specification

Including:

- High integrity design - less than 0.05% vol/hr leak rate
- Modular construction
- Positive or negative pressure
- System alarms to requirements
- Specifications and options available to suit customer's individual needs



### Optional features

Options include:

- ULPA filter
- Electrical sockets
- Service taps
- Bag out ports
- Carbon column
- Modular side extension of glovebox system
- Sleeve test kit

Further options are available on request

### Other products and services

Bassaire has over fifty years experience in the design, manufacture, installation and maintenance of cleanrooms and clean air products.



cleanrooms



laboratories



clean air products



servicing

ds-isolator-v2

[www.bassaire.co.uk](http://www.bassaire.co.uk)

www.bassaire.co.uk | Tel: +44 (0) 1489 885111 | Fax: +44 (0) 1489 885211 | Email: sales@bassaire.co.uk

Bassaire Limited · Duncan Road · Park Gate · Southampton · Hampshire · SO31 1ZS · England

Registered in England: 2107555 · © Bassaire Limited